

Alu Lab Student Guide

Introduction to PCR — The Polymerase Chain Reaction

You are about to perform a procedure known as PCR¹ to amplify a specific sequence of your own DNA in a test tube. You will be looking for a particular piece of DNA that is present in the genes of many, but not all, people. Analysis of the data generated in this laboratory will enable you to determine whether or not you carry this specific DNA sequence.

The **genome**, composed of DNA, is our hereditary code. This is the so-called blueprint that controls much of our appearance, behavior, and tendencies. **Molecular biology** is the study of genes and the molecular details that regulate the flow of genetic information from DNA to RNA to proteins, from generation to generation. **Biotechnology** uses this knowledge to manipulate organisms' (microbes, plants, or animals) DNA to help solve human problems.

Within the molecular framework of biology, DNA, RNA, and proteins are closely tied to each other. Because proteins and enzymes ultimately play such a critical role in the life process, scientists have spent many lifetimes studying proteins in an attempt to understand how they work. With this understanding, it was believed we could cure, prevent, and overcome disease and physical handicaps as well as explain exactly how and why organisms exist, propagate, and die. However, the complete answer to how and why does not lie solely in the knowledge of how enzymes function; we must learn how they are made. If each enzyme is different, then what controls these differences and what is the blueprint for this difference? That answer lies within our genome, or genetic code.

Thus, you may realize why researchers today, in an attempt to understand the mechanisms behind the various biological processes, study nucleic acids as well as proteins to get a complete picture. In the last 20 years, many advances in nucleic acid techniques have allowed researchers to study the roles that nucleic acids play in biology. It took the imagination and hard work of many scientists to reveal the answers to one of the most mysterious puzzles of life — understanding the mechanisms that control how DNA is translated into proteins within living cells.

Before beginning this lab, think about your answers to the following questions

How is DNA faithfully passed on from generation to generation? What causes genetic disease in some people but not others? How do scientists obtain DNA to study? What secrets can DNA tell us about our origins? What human problems can an understanding of DNA help us solve? Should we unlock the secrets held in this most basic building block of life?

PCR Set the Stage for a Scientific Revolution

In 1983, Kary Mullis² at Cetus Corporation developed the molecular biology technique that has since revolutionized genetic research. This technique, called the **polymerase chain reaction** (PCR), transformed molecular biology into a multidisciplinary research field within 5 years of its invention. Before PCR, the molecular biology techniques used to study DNA required such a high level of expertise that relatively few scientists could use them.

The objective of PCR is to produce a large amount of DNA in a test tube (in vitro), starting from only a trace amount. Technically speaking, this means the controlled enzymatic amplification of a DNA sequence, or gene, of interest. The template strands can be any form of DNA, such as genomic DNA. A researcher can use tiny amounts of genomic DNA

from a drop of blood, a single hair follicle, or a cheek cell, and make enough DNA to study. In theory, only a single template strand is needed to copy and generate millions of new identical DNA molecules. Prior to PCR, this would have been impossible. It is the ability to amplify the precise sequence of DNA of interest that is the true power of PCR.

PCR has made an impact on four main areas of genetic research: gene mapping; cloning; DNA sequencing; and gene detection. PCR is now used as a medical diagnostic tool to detect specific mutations that may cause genetic disease;³ in criminal investigations and courts of law to identify suspects,⁴ and in the sequencing of the human genome.⁵ Prior to PCR, the use of molecular biology techniques for therapeutic, forensic, pharmaceutical, agricultural, or medical diagnostic purposes was neither practical nor cost-effective. The development of PCR transformed molecular biology from a difficult science to one of the most accessible and widely used disciplines of **biotechnology**.

Lesson 1 Cheek Cell DNA Template Preparation

To obtain DNA for use in the polymerase chain reaction (PCR) you will extract the DNA from your own living cells. It is interesting to note that DNA can be also extracted from mummies and fossilized dinosaur bones. In this lab activity, you will isolate DNA from epithelial cells that line the inside of your cheek. To do this, you will rinse your mouth with a saline (salt) solution, and collect the cells using a centrifuge. You will then boil the cells to rupture them and release the DNA they contain. To obtain pure DNA for PCR, you will use the following procedure:

The cheek cells are transferred to a micro test tube containing **InstaGene™ matrix**. This particulate matrix is made up of negatively charged, microscopic beads that chelate, or grab, metal ions out of solution. It traps metal ions, such as Mg^{2+} , which are required as catalysts or **cofactors** in enzymatic reactions. Your cheek cells will then be **lysed**, or ruptured, by heating to release all of their cellular constituents, including enzymes that were once contained in the cheek-cell lysosomes. Lysosomes are sacs in the cytoplasm that contain powerful enzymes, such as **DNases**, which are used by cells to digest the DNA of invading viruses. When you rupture the cells, these DNases can digest the released DNA. However, when the cells are lysed in the presence of the chelating beads, the cofactors are absorbed and are not available to the enzymes. This virtually blocks enzymatic degradation of the extracted DNA so you can use it as the template in your PCR reaction.

You will first suspend your isolated cheek cells in the InstaGene matrix and incubate them at 56°C for 10 minutes. This preincubation step helps to soften plasma membranes and release clumps of cells from each other. The heat also inactivates enzymes, such as DNases, which can degrade the DNA template. After this 10 minute incubation period, place the cells in a 100°C water bath/heat block for 5 minutes. This high temperature ruptures the cells and releases DNA from their nuclei. You will use the extracted genomic DNA as the target template for PCR amplification.

Lesson 1 Cheek Cell DNA Template Preparation (Lab Protocol)

Workstation Checklist

Materials and supplies required at the workstations prior to beginning this exercise are listed below.

Student Workstation	Quantity per Station	(✓)
1.5 ml microfuge tubes	4	<input type="checkbox"/>
Permanent markers (variety of colors)	2	<input type="checkbox"/>
Cups with 10 ml 0.9% saline	4	<input type="checkbox"/>
Waste container containing bleach solution	1	<input type="checkbox"/>
P-1000 micropipette, 200–1000 μ L	2	<input type="checkbox"/>
Pipette tips, blue 200–1000 μ L	8	<input type="checkbox"/>
Screw cap tubes with 200 μ L InstaGene matrix	4	<input type="checkbox"/>
Copy of Lab Procedure	4	<input type="checkbox"/>
Epitube rack	1	<input type="checkbox"/>


Class Equipment	Quantity per Class	
Heat block (56°C)	1	<input type="checkbox"/>
Heat Block (100°C)	1	<input type="checkbox"/>
Centrifuge	4	<input type="checkbox"/>
Vortexer (optional)	4	<input type="checkbox"/>

Lesson 1 Cheek Cell DNA Template Preparation (Lab Protocol)

- Each member of your team should have a cup containing 10 ml of 0.9% saline solution. Label your cup with a 4 digit pin number that only you know.



- (Note: Do not throw away the saline until you complete step 6 below.) Pour the saline from the cup into your mouth. Rinse vigorously for 30 seconds. While rinsing, grip the inside of your cheeks with your teeth to help free old cells. Expel the saline back into the cup without anyone seeing you do this step. –Don't look into anyone else's cup.




- Set a P-1000 micropipette to 1,000 μL and transfer 1 ml of your oral rinse into the microfuge tube labeled with your 4 digit pin number. Use the permanent marker to add colorful identifying marks to the lid and hinge of your microfuge tube.



4.




- Spin your tube in a balanced centrifuge for 2 minutes at full speed (10,000 x g). When the centrifuge has completely stopped, remove your tube. You should be able to see a pellet of whitish cells at one side near the bottom of the tube. Ideally, the pellet should be about the size of a match head. Compare your pellet size to your lab group members. If you can't see your pellet, or your pellet is too small, pour off the saline supernatant, add more of your saline rinse from your cup, and spin again.



2 minutes


- Pour the supernatant (liquid above a solid) back into your disposable cup. Take care not to lose your cell pellet, carefully blot your micro test tube on a tissue or paper towel. It's ok for a small amount of saline to remain in the bottom of the tube (as you see in the drawing at the right).




- Re-suspend the pellet thoroughly by vortexing, flicking, or racking the tube until no cell clumps remain.




7. Obtain a screw cap tube containing 200 μ l of InstaGene. Compare the amount of InstaGene beads (bottom of tube) to your lab partners. (The InstaGene beads are heavy and tend to fall to the bottom of a solution which makes dividing the source solution challenging.) Label a screw cap tube containing 200 μ L of InstaGene with your 4 digit pin number and colorful markings on and around the lid and sides. Adjust the P-1000 Micropipette to 200 μ l. Using the pipette transfer all of your re-suspended cells into the screw cap tube containing InstaGene. You may need to use the pipette a few times to transfer all of your cells.
8. Screw the cap tightly on the tube. Shake, rack, or vortex to mix the contents.
9. Place your screw cap tube into the 56°C heat block for 10 minutes. At the halfway point (5 minutes), shake or vortex the tubes several times. Place the tubes back in the heat block for the remaining 5 minutes.



10. Remove the tubes from the heat block and shake them several times. Now place your tube in a slot of the heat block at 100°C for 5 minutes.



11. Remove the tubes from the 100°C heat block and shake, rack or vortex several times to re-suspend the sample. Place your screw cap tube in a balanced arrangement in a centrifuge (because our screw cap tubes vary so greatly in volume of liquid and InstaGene, only the old centrifuge works for this step). Pellet the matrix by spinning for 5 minutes at 6,000 x g (or 10 minutes at 2,000 x g).



12. Store your screw cap tube with the rest of the class in a rack in the refrigerator until the next laboratory period.

CLEAN UP: Pour the contents of your paper cup directly down the drain of the sink and throw your paper cup in the garbage. Open your empty flip top microfuge tube and place it in the bleach solution. Organize your lab station including stools under the table.

BACKGROUND Lesson 1 Cheek Cell DNA Template Preparation

Processing Cheek Cell Samples to Obtain a Genomic DNA Template for PCR

A. InstaGene matrix: What function does it perform?

InstaGene matrix consists of a suspension of negatively charged Chelex® microscopic beads, which bind divalent cations like magnesium (Mg_{2+}). It is important to remove divalent cations from genomic DNA samples because the cations assist enzymes that degrade the DNA template. When cheek cells or hair follicles are lysed and boiled in the presence of InstaGene matrix, the divalent cations released from the cells bind to the beads, and the heat inactivates the DNA-degrading enzymes. The beads are then pelleted by centrifugation. The supernatant (light liquid above the heavier tube contents), which contains clean, intact genomic DNA, can be used as the template in the PCR reactions.

Each student will prepare genomic DNA from cheek cells isolated using a saline mouthwash or from hair follicles. For students using the cheek cell protocol, 1 ml of cells collected using the mouthwash should provide sufficient material for DNA preparation. Some students may need to use 2 ml or more of the saline mouthwash to obtain sufficient cells to prepare DNA. Please note: it is **not** recommended to use more than 3 ml of the saline mouthwash to prepare DNA, (see Troubleshooting below). Once the cells have been spun in a centrifuge, a cell pellet about the size of a match head should yield enough cells for subsequent steps. Eating just prior to cell collection is not recommended, as food particles may make cell preparation more difficult.

Troubleshooting

Multiple explanations can account for poor student results

1. **Inadequate collection of cheek cells.** A visible cell pellet about the size of a match head should be obtained following centrifugation of the saline mouthwash. If no cell pellet is visible, or the pellet is too small, additional swished saline may be centrifuged until a pellet of the desired size is obtained. However, it is not recommended to collect more than 3 ml of cells (see point below).
2. **Excessive number of cells.** Just as too few cells will yield insufficient genomic DNA, an excessive number of cells will saturate the capacity of the InstaGene, resulting in little or no amplification.
3. **InstaGene matrix not transferred.** Each workstation is supplied with tubes of InstaGene matrix that were aliquotted by the instructor and placed on ice. These tubes of matrix must be mixed prior to each pipetting to bring the beads up into suspension. If no beads were transferred into the student's tube, the divalent cations will not be removed from the genomic DNA preparation, and the PCR reaction will be inhibited.

If the DNA samples will not be amplified within 24 hours, they can be stored in the refrigerator in the InstaGene matrix for up to 1 week. For longer storage, place samples in the freezer to prevent DNA degradation. Before the samples are used in PCR, the beads should be re-pelleted by centrifugation just prior to making up the PCR reactions. However, processing the samples within 24 hours is recommended.

What are the functions of each incubation step?

The preincubation step is carried out at 56°C and performs two functions:

1. Heating the cell suspension aids in the breakup of connective tissue that holds the cells together. Breaking up the tissue makes the cells easier to lyse during the subsequent 100°C incubation step.
2. Preincubation at 56°C inactivates DNases, enzymes that are naturally present in the cell suspensions, and which could degrade the genomic DNA and inhibit the PCR reactions.

Heating the cell samples to 100°C ruptures the cell membranes, thereby releasing the cellular contents, which include genomic DNA. The genomic DNA serves as the template in the PCR reactions.

Chapter 12 Alu Lab

Lesson 1 DNA Template Preparation

Focus Questions

Name _____ Per__

1. How would you make approximately 500mL of a 0.9% salt solution? Show your work
2. Why do you think we are using a salt solution for the mouthwash step?
3. Thinking back to your biotech experience in Biology, write down easy to follow instructions for how to use a micropipette. (Include how to withdraw a sample from a microfuge tube (reaction tube) taking it up into the micropipette, and how to expel the sample from the micropipette into a reaction tube.
4. According to the lab, which micropipettes will we be using in step 3 and 7? What is the range of each (look back in your notes that we took about the micropipette)? From top down, which three digits will be displayed on the micropipette in step 3 and 7?
5. Why do you think you do not pour out the saline right after step 3?
6. At step 4 what is the purpose of using the centrifuge for 2 minutes at high speed?
7. At step 5, what is in the supernatant? What is in the pellet?
8. What do steps 5 and 6 accomplish?
9. What is the purpose of the 56°C hot water bath and the 100°C heat block?
10. Why do you think we are using screw-cap tubes instead of the usually flip top tubes for these steps?
11. Why is it necessary to chelate the metal ions (using InstaGene matrix) from solution during the boiling/lysis step at 100°C? What would happen if you did not use a chelating agent such as the InstaGene matrix?
12. What is needed from your cells for PCR?
13. What do you think is the purpose of the centrifuge at step 11?
14. Why do you think the DNA is stored cold with the InstaGene matrix after boiling the samples?
15. Be sure you are able to describe why we are doing each step of the lab procedure.

Informed Consent

Date: January 2019

Dear Parent(s) or Guardian,

The AP Biology class at Lynbrook High school has the opportunity to participate in a class exercise in which an important technique in biotechnology will be used to analyze the students' DNA. The technique the students will be using is called the Polymerase Chain Reaction (PCR). It is a method by which a particular piece of DNA can be amplified many million-fold. PCR has a number of applications in the scientific community, including uses in forensics, diagnostics, parentage testing, and evolutionary studies.

In this laboratory protocol, students will be isolating DNA from their own cheek cells. They will then apply the PCR technique to amplify a particular segment of their DNA. This segment is not known to be associated with any genetic disease and variation between individuals in this region is in no way an indicator of health or genetic fitness. The results of this particular lab exercise are for teaching purposes only and will NOT be used for any diagnostic or identification purposes.

Your student's privacy will be protected. The student's name will not be linked to his/her DNA.

Participation is voluntary. By signing this permission form, you are allowing your student to use their own DNA sample in this lab investigation. If you would prefer that your student not use their own DNA sample, your student will still participate in each step of the lab procedure by using one of the provided DNA controls.

Sincerely,

Ms. Della Santina

nicole_dellasantina@fuhisd.org

I _____ give permission for my student _____ to
print parent/guardian name print student name

use their own DNA sample for the AP Biology lab.

Parent/Guardian Signature of Consent

Date

Informed Consent

Date: January 2019

Dear Parent(s) or Guardian,

The AP Biology class at Lynbrook High school has the opportunity to participate in a class exercise in which an important technique in biotechnology will be used to analyze the students' DNA. The technique the students will be using is called the Polymerase Chain Reaction (PCR). It is a method by which a particular piece of DNA can be amplified many million-fold. PCR has a number of applications in the scientific community, including uses in forensics, diagnostics, parentage testing, and evolutionary studies.

In this laboratory protocol, students will be isolating DNA from their own cheek cells. They will then apply the PCR technique to amplify a particular segment of their DNA. This segment is not known to be associated with any genetic disease and variation between individuals in this region is in no way an indicator of health or genetic fitness. The results of this particular lab exercise are for teaching purposes only and will NOT be used for any diagnostic or identification purposes.

Your student's privacy will be protected. The student's name will not be linked to his/her DNA.

Participation is voluntary. By signing this permission form, you are allowing your student to use their own DNA sample in this lab investigation. If you would prefer that your student not use their own DNA sample, your student will still participate in each step of the lab procedure by using one of the provided DNA controls.

Sincerely,

Ms. Della Santina

nicole_dellasantina@fuhisd.org

I _____ give permission for my student _____ to
print parent/guardian name print student name

use their own DNA sample for the AP Biology lab.

Parent/Guardian Signature of Consent

Date